Appareils de laboratoire

Chauffage

Laboratory equipment

Heating

Réf :
701 311

Français – p 1

English – p 11

Version : 4104

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Sérénit[®] laboratory electric burner

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1	Description	8
2	Characteristics	8
3	Warning symbol	8
4	Starting up	9
4.1	First use	9
4.2	Starting up	9
5	Safety	9
5.1	Thermal safety	9
5.2	Electrical safety	10
6	Precautions for use	10
6.1	Handling the device	10
6.2	Protection of persons	10
6.3	Gas fumes	10
6.4	Precautions for use	11
6.5	Working area	11
6.6	Overturning of the device	11
6.7	Warning	11
7	Use	11
7.1	Manual heating of test tubes	11
7.2	Performing static fittings:	12
7.3	Combustion and production of flames:	13
8	Storage	13
9	Care and maintenance	14
10	Accessories	14
10.1	Accessories supplied	14
10.2	Optional accessories	14
11	Warranty and after sales service	14

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1 Description

The Sérénit® laboratory electric burner is a heating device that can replace the Bunsen burner to heat test tubes, containers and to perform various chemistry experiments that require the use of a flame.


Diagram 1

2 Characteristics

- Mass: 2.6 kg.
- Dimensions (mm): **L.280 x w.180 x h.125.**
- Removable rod stand, length: 300 mm.
- Power: 500 W
- Thermostatic setting of the temperature by energy regulator.
- Protection against overloads by a fuse (type T 3,15 A, L 250 V).
- Protection against abnormal overheating by non-resettable thermal protector calibrated at 100 °C.
- Standardised cord of 1.90 m length
- Power supply: 230 V (2 P+T) 50/60 Hz.
- Category of use II.

Environmental conditions:

- Operating temperature between 15 and 35 °C.
- Relative humidity less than 75 %.

3 Warning symbol

This symbol means « CAUTION » and Asks the user to thesis Refer to instructions.

This symbol means « CAUTION » hot surface.

4 Starting up

4.1 First use

When operating for the first time or following storage for a long duration in an atmosphere with a high humidity level, the Sérénit® electric burner can produce leakage currents and very rarely cause the earth leakage circuit breaker to trip. This problem is caused by the very nature of the heating resistance used in this type of device and **does not present any danger to the user** as long as the device is connected to an electrical socket with a grounding terminal (Class I device). To solve this problem, power on the devices in position 10 for about 45 min in order to recondition the heating resistance. To perform this operation, it is important not to connect more than two devices to each 30 mA earth leakage circuit breaker.

4.2 Starting up

The device will be connected to a mains supply corresponding to its electrical characteristics, with earthing and differential protection.

After connecting to the mains, the power on indicator lamp (1) comes on (red indicator lamp).

Position the setting mark of the button (3) on the desired temperature; the heating indicator lamp (2) comes on when the heating resistance is powered by the electronic dimmer switch of the device.


Diagram 2

5 Safety

5.1 Thermal safety

The Sérénit® electric burner is designed to provide the user with maximum protection against any risk of accidental burns.

The heating area is delimited by a double ring and a safety cradle to prevent any accidental contact with hot parts. The safety cradle with a deliberately inclined shape does not allow large containers to be placed on the device as these would obstruct lateral ventilation (Diagram 3).

Diagram 3


Diagram 4

The Sérénit® electric burner is designed based on an open structure ensuring passive thermal insulation so that the lateral walls of the device remain close to the ambient temperature even after a long period of use (diagram 4).

This open structure additionally helps the flow of spilt liquids through the device for increased safety.

Diagram 5

Diagram 6

A removable grille system (Part no. 703 530) is used to effectively complement this protection device for greater operating safety (diagrams 5 and 6). This grille also prevents the containers from toppling.

5.2 Electrical safety

The Sérénit® electric burner is equipped with a shielded element, a protection fuse and a non-resettable thermal circuit breaker that can be used to switch off the circuit in case of abnormal heating of the device.

This device complies with the standards:

- EN 61010-1 ed 2010 et EN 61010-2-010 ed 2003.

The device is powered off by removing the mains plug, which must remain accessible.

6 Precautions for use

6.1 Handling the device

The Sérénit® laboratory electric burner is a heating device that must be handled carefully.

It is recommended not to move the device when it is heating or powered on.

Hold the device only by the non-metallic parts whose temperature is close to the ambient temperature.

6.2 Protection of persons

Never carry out any experiment without acquainting yourself beforehand with the risks relating to the chemicals used.

Work with personal and collective protective equipment adapted to the risks inherent in the chemicals used (Hood, Goggles, etc.).

Do not use this device to heat flammable liquids, products or product mixtures that present the risk of ignition or explosive decomposition.

6.3 Gas fumes

Work in a ventilated area if there is a risk of gas fumes during experiments performed.

6.4 Precautions for use

Never let the device operate without supervision.
Never leave an empty container on the device when heating.
Use only borosilicate glass containers to perform experiments (metal containers are prohibited).

6.5 Working area

Never use the laboratory electric burner on a flammable surface.
Never work at less than 30 cm from a wall.

6.6 Overturning of the device

In order to prevent the device from overturning, it is recommended to:

- Fasten heavy attachments on the removable rod stand
- Hold the device by the removable rod stand.

Never replace the removable rod stand supplied with the device by another rod.

6.7 Warning

The JEULIN Company warns the user that the protection of this device may be compromised if it is used in a manner that is not specified in these instructions.

7 Use

7.1 Manual heating of test tubes

Diagram 7

Diagram 8

It is possible to heat the test tubes by placing them above the heating area to control the reaction accurately. The risks of the spattering of liquids when boiling may be avoided by heating the tube over its entire length by tilting it over the heating area (diagram 7).

To perform experiments requiring a strong heat input such as pyrolyses or production of metal alloys, it is possible to insert the test tube (up to 20 mm diameter) into the opening located at the centre of the grille so that the tube is as close to the heating resistance as possible (diagram 8).

7.2 Performing static fittings

The Sérénit® electric burner is designed to heat containers and enable fittings.

To do this, the Sérénit® electric burner is fitted with a removable rod stand that can be placed at the 4 corners of the device by screwing it to the angle pads (diagram 9).

Diagram 9

Diagram 10

Diagram 11

Diagram 12

It is thus possible to safely heat containers of up to 105 mm diameter such as beakers, flasks and conical flasks. The safety cradle in addition ensures the centring of the containers over the heating area (diagrams 10, 11 and 12).

Diagram 13

Diagram 14

The removable rod stand of the Sérénit® electric burner is used to carry out complex assemblies such as distillation or temperature measurement set-ups (Diagram 13 and 14).

7.3 Combustion and production of flames

The heating area is designed to enable combustion of various material such as mineral charcoal, steel wool, magnesium, etc. (Diagrams 15, 16 and 17).

Diagram 15

Diagram 16

Diagram 17

The ignition is initiated by mere contact with the heating resistance of the device.

It is also possible to safely produce a flame on the heating area of the Sérénit® electric burner (Diagram 18) by using special combustible pellets (Part no. 101 012).

Diagram 18

Diagram 19

Producing a flame is very useful for combustion of metal powders and flame colour tests (Diagram 19).

8 Storage

To make handling easier, the device is equipped with a grasping handle located under the control console.

The shape of the devices (diagram 20) is designed to enable them to be piled up (This operation is possible only when the removable rod stand is removed).

Diagram 20

9 Care and maintenance

This device does not require periodic or preventive maintenance. In case of pollution, this device will have to be cleaned according to the following precautions:

- **The device must be unplugged from the mains supply.**
- Any cleaning (whatever be the product used) must be done only after the device has cooled down.

Check that the lower discharge hole of the bowl is not blocked.

Any incinerated matter will have to be removed from the grille and the resistance.

This device does not require any maintenance of the electrical part.

In case of failure:

If the red LED Power does not turn on:

- Set the knob to O
- Disconnect the mains
- Change the fuse on the side of the device: fuse T3, 15 A, L 250 V.

If, after checking the fuse, the unit still does not work, please contact our company.

This device must in no case be dismantled by the user for any reason whatsoever. Any intervention will have to be done by the JEULIN Company or its representative.

10 Accessories

10.1 Accessories supplied

The device is delivered with a removable rod stand.

This rod can be placed on the four angle pads located around the heating area. It is fixed simply by screwing.

10.2 Optional accessories

Removable protective screen
Part no. 703 530

Polynux® Nut
Part no. 703 452

Combustible pellets
(box of 60)
Part no. 101 012

11 Warranty and after sales service

The warranty is 2 years, the equipment must be returned to our workshops.

For any repairs, adjustments or parts, please contact us.

All manipulations of the device by an unauthorized person automatically cancels the warranty period and may compromise the safety of the user.

JEULIN – TECHNICAL SUPPORT
468 rue Jacques Monod
CS 21900
27019 EVREUX CEDEX FRANCE
+33 (0)2 32 29 40 50

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Direct connection for technical support

A team of experts
at your disposal
from Monday to Friday
(opening hours)

- You're looking for technical information ?
- You wish advice for use ?
- You need an urgent diagnosis ?

We take in charge your request
immediatly to provide you
with the right answers regarding
your activity field : Biology, Physics,
Chemistry, Technology.

Free service*

+33 2 32 29 40 50**

* Call cost not included.
** Only for call from foreign countries.

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SAS au capital de 1 000 000 € - TVA intracommunautaire FR47 344 652 490 - Siren 344 652 490 RCS Evreux