

Technical Datasheet

Article No. 2941

Peroxide 1000 - Test Strips

0 – 1000 mg/L, 100 strips, semi-quantitative

Single-use test strips for the rapid, semi-quantitative detection of peroxide residues in solutions.

Parameter	Value
Measuring range	0 - 1000 mg/l
Gradation	0 - 100 - 200 - 400 - 600 - 800 - 1000 mg/l
Colour reaction	White to dark brown
Shelf life	3 years from the date of manufacture Use opened packages within 3 months

Instructions for Use

1. Dip the test strip into the test solution for 2 seconds.
2. Shake off any excess solution.
3. Wait 30 seconds.
4. Compare the color that appears on the test strip with the color chart printed on the package.

Storage Instructions

Unused sticks must be stored in their original packaging. The packaging must be kept tightly sealed in a cool, dry place where it is not exposed to chemical fumes.

Recommended storage temperature: 2 – 8 °C, storage at higher temperatures > 8 °C might shorten the shelf-life.

Version No. 01/2026-05-11